



Features

- FAST CHAIR PREMIUM POWERTRAXX w/PLASTIC SEAT (Ordering No. 60-4170-001) for use as a stair, evacuation, or transport chair
- Cost-effective solution featuring Ferno Powertraxx track system
- Stair track easily attachable or removable to minimize operator carrying weight
- Lightweight all-aluminium frame supports patients up to 250 kg
- Locking in the folded position
- Comfortable locking and unlocking
- Anti-tip protection for the main frame
- Side handle protective guards
- Smooth and quiet wheels - two 125 mm front swivel and two 150 mm rear fixed with locks
- Wheels have minimum debris collection while maintaining good traction and comfort, and provide very good maneuverability in confined spaces such as stair landings or tight corridors
- Head support, three restraints, comfortable armrests and comfort footbar included
- Wide rear locking handles
- Telescoping handle is extendable and adjustable in four positions (1000 to 1580 mm)
- Front telescoping handles adjustable in five positions (maximum extension of 340 mm)
- Plastic cover is easy to clean

Specifications

Height min	Overall Width	Overall Length	Seat Width	Folded L,W, Thickness	Weight with track and battery	Load Limit
1020 mm	575 mm	725 mm	420 mm	1020x575x215 mm	30 kg	250 kg

Accessories

FAST CHAIR WALL MOUNT

designed to securely hold the FAST Chair in both a static or dynamic application. Single hand release allows for quick chair fastening or detachment.

60-4190-001



4- POINT RESTRAINT for enhanced safety

62-4400-001



FAST BIG WHEELS

- suitable for off-road use
- can be quickly and easily attached

65-0407-001-P



MESH BASKET

for storing accessories such as oxygen bottle

65-4123-001



The FAST CHAIR PREMIUM POWERTRAXX with plastic seat complies with and exceeds the requirements of the international standard EN 1865-6:2024.

