

UNIVERSAL Stair Chair **FASTENER**



0567086 - Fastener without ICS

0567087 - Fastener with ICS

The **FERNO® Universal Stair Chair Fastener** is designed to secure your compatible stair chair to the wall of your ambulance. Features an adjustable lock handle and wheel blocks.




Integrated Charging System*

Charge your chair during transport or anytime.



Lock handle keeps your chair secure. Quick, one-handed release.



Wheel blocks for improved stability and retention

*Designed to work with the EnduraCharge Battery



UNIVERSAL

Stair Chair

FASTENER

Specifications

HEIGHT	37.5 in 950 mm
WIDTH	24 in 607 mm
DEPTH	11.75 in 300 mm
WEIGHT (Fastener Only)	22 lb 10 kg
WEIGHT (Fastener with ICS)	27 lb 12.2 kg



Stair Chair not included

The Universal Stair Chair Fastener is compatible with the following chairs:

- Ferno Transcend series (excluding the Evac Chair)
- Ferno EZ Glide series
- Ferno Venice series
- Ferno FAST Chair series
- Stryker Stair Pro