



**Wheel-Kit 71 W-V for
FERNO Basket Stretcher 71/71 S**

WHEEL-KIT 71 W-V



Wheel-Kit 71 W-V for FERNO Basket Stretcher 71/71 S

The wheel kit for our basket stretcher 71 or 71 S facilitates transport over longer distances.

With a few simple steps, the kit is attached to the foot part of the basket stretcher with 4 quick-release fasteners.

The 250 mm big solid rubber wheels ensure easy running even on uneven terrain.

Powder coated aluminum frame.



GENERAL SPECIFICATIONS

	Length	Width	Height	Weight	Load Capacity
Wheel-Kit 71 W-V	600 mm	680 mm	320 mm	5,1 kg	280 kg