



SCOOP EXL[®] TSL EXTENDER KIT

The Ferno[®] Scoop[®] Extender kit TSL is an innovative accessory designed to expand the patient surface on your Scoop[®] EXL. This extension allows the width of the Scoop[®] EXL to be adjusted to match the patient's physical size, increasing flexibility and ease of use.

With a safe load capacity of 227 kg and adjustable length, Ferno's Scoop[®] EXL can now be used on patients of all ages and sizes. You can choose to expand the width only at the head end for patients with broad shoulders or a strong upper body, or extend both the head and foot ends for a generally wider lying surface.

The extender kit is designed to be quick and easy to use. The TSL clicks securely into place in the lock, allowing the Scoop[®] EXL to be adjusted in seconds without hassle.



KEY FEATURES

- Easy to take with you
- Flexible use
- Can be used in one or both locks, all depending on the situation and patient size
- One set contains 2 pcs. TSL to the locks
- Load capacity still 227 kg
- Provides a better fit if the Scoop stretcher is combined with Lucas for cardiac arrest

Art. Nr	F0832274
Certifications	MDR
Load Capacity	227 kg
Material	Metal
Weight	490 g
Height	35 mm
Length	130 mm
Width	35 mm

