



Features

- Independent leg operation provides safe loading into the ambulance, legs unfold automatically for easy unloading
 - Color-coded controls are easy to use and require minimal training
 - F3 Monobloc is easy to handle, supports up to 250kg and has a strong structure in aluminum
 - Four 200mm diameter swivel wheels with brakes
 - Two front wheels with swivel lock
 - Features raised or contour shock position
 - Swing-down sidearms
 - Four telescopic carrying handles with five positions (longer and shorter included)
 - Four 100mm diameter wheels for the stretcher
 - Foot Support
 - Secure hardboard (BB) insertion system placed under the cot
 - The rugged and lightweight IV Pole offers an integrated height adjustment which allows operators to set their desired height for treatment
 - F3 has been designed to comply with and exceed international Standards such as EN 1865-1:2010+A1:2015 and EN 1789:2020
-
- 60-0340-001 MONOBLOC F3 4P 59 cm (H56)
 - 60-0340-002 MONOBLOC F3 4P 54 cm (H50)
 - 60-0340-003 MONOBLOC F3 4P 69 cm (H69)

Accessories & Fasteners



F2 LOCKING DEVICE



MATTRESS
ORANGE/BLACK



RESTRAINTS KIT



HEADREST



INSTRUMENT TRAY

Specifications

Length min.	Length max.	Width	Load Height 1	Load Height 2	Weight	Load limit
1950 mm	2350 mm	568 mm	500 mm	540 mm	45 kg	250 kg
			560 mm	590 mm		
			650 mm	690 mm		