

Patient Cover - Basket Stretcher fitted Elastic with Zipper

Design

The Traverse Rescue Patient Cover protects the patient from rain, wind, sleet and snow while still giving easy access through the double zipper centre opening to monitor the patient's condition or treat their injuries. Made from coated nylon and held in place by an elastic rim and Velcro attachment straps, the Patient Cover can be rolled up and stuffed into it's own storage pouch which can be attached to the stretcher with the Velcro attachment strap. Colour: Royal Blue



Part No.	Description
15-0130	Patient Cover - Basket Stretcher fitted Elastic with Zipper